

Beyes®

The Next Standard



November
26th-29th, 2023
BOOTH #4400
SPECIALS

AirLight® M800 Plus

LED High Speed Handpiece

BUY 3 FOR \$1100!

- 23W of Torque
- 410,000 RPM
- Ceramic Ball Bearings
- Handles up to 60 PSI of air pressure
- Titanium Body
- 1 Year Warranty

Backend	Code
4 Hole	HP3008P
Beyes®	HP6018P
W&H®	HP3022P
Star®	HP3054P
NSK®	HP3043P
KAVO®	HP3036P



Reg. \$548/ea



Maxso®

Low Speed Package

- 1 Year Warranty

Includes:

Contra Angle-Push Button (SL2026)
Straight Nose Cone (ST2010)
Air Motor-NS/M4 Backend (MT2010)



Reg. \$470/pk

Sale: \$439 /pk

Maxso® E800

Electric Handpiece System

- **Powerful:** 5.0 Ncm of torque
- **Steady Performance:**
Speed of rotations is controlled between 2,000 and 40,000min-1
- **Smooth Operation:**
Torque ripple and heat reduction technology eliminates vibrations
- **Convenient & Durable:**
Smooth and easy to clean surface with autoclavable brushless motor
- **3 Programmable Memory Settings:**
Easily switch between speed presets for more efficient workflow
- **New Touch Screen Interface:**
Easily and quickly navigate through menu settings on a color LCD display panel

Code: MT8001 Reg. \$1,678

MT8001KIT1: E800P + One X99L Attachment

Reg. \$2,468 **Sale: \$2,075**

MT8001KIT2: E800P + Two X99L Attachment

Reg. \$3,202 **Sale: \$2,699**



Attachment:

BUY 2 GET 1 FREE!

Maxso® X99L

1:5 Contra Angle

- Anti-Heat System
- Ceramic Ball Bearing
- ISO E-Type Connection
- Quattro Spray
- Fiber-optic
- 1 Year Warranty

Code: SL3025

Reg. \$992/ea **Net: \$661.33/ea**



Maxso® M200E

Non-Optic High Speed Handpiece

BUY 5 FOR \$550!

- 19W of Torque
- 400,000 RPM
- Single Spray
- Push Button
- Low Turbine Noise
- No Swivel
- 6 Months Warranty



Backend	Code
4 Hole	HP2105

Reg. \$128/ea

Duray® Intraoral Sensor Package

Includes size 1 & size 2 kits

- Thinnest Sensor 4.4mm
- Durable and Reliable
- High Resolution Imaging
- Easy To Use
- 24 Months Warranty

Code: DS3001

Reg. \$6,298 **Sale: \$5,250**

Size 2: DS2001

Reg. \$3,676 **Sale: \$3,125**

Size 1: DS1001

Reg. \$2,940 **Sale: \$2,450**



Slimax-C Plus®

LED Curing Light System

BUY 1, GET 1 FREE!

- 5-Second Curing
- 2000 mW/cm ± 10%
- Built-In Digital Radiometer
- Cordless Design
- Rechargeable Samsung Battery
- 1 Year Warranty



Color	Code
Silver	SM1003P-S
Blue	SM1003P-B
Pink	SM1003P-P

Reg. \$660/ea
Net: \$330/ea



Duray® Imaging Plate Scanner Package

DC1, Duray Scanner DC1 Package

- Powerful Processing
- Smart Design
- Compact Tabletop Design
- Comfortable
- 24 Months Warranty

Code: DC1001

Reg. \$6,826 **Sale: \$5,550**



Airlight® M800X

LED High Speed Handpiece

- 27W of Torque
- 430,000 RPM
- Instant-Stop Capability
- X-Ball Bearings
- Extra Quiet
- Micro-Tex Slip-Proof Coating
- 25 Months Warranty

Backend	Code
Beyes®	HP6018X
Star®	HP3054X
NSK®	HP3043X
KAVO®	HP3036X

Reg. \$1,275/ea **Net: \$637.50/ea**

BUY 1, GET 1 FREE!



Airlight® M800 Surgical 45

LED High Speed Handpiece

- 22W of Torque
- 45 Degree Angled Head
- 410,000 RPM
- Rear Air Exhaust
- Direct LED
- Stream Water
- 1 Year Warranty

Backend	Code
4 Hole	HP3048
Beyes®	HP6038
NSK®	HP3044
KAVO®	HP3037

Reg. \$598- **Sale: \$498/ea**



Maxso® X200X

Fiber-Optic High Speed Handpiece

- 27W of Torque
- 430,000 RPM
- Instant Stop Capability
- Quattro Spray
- Titanium Body
- 25 Month Warranty

Backend	Code
NSK®	HP2066X
KAVO®	HP2060X

Reg. \$1,048/ea **Net: \$524/ea**

BUY 1, GET 1 FREE!



Maxso® X92

20:1 Reduction Contra Angle

- 20:1 Ratio
- Non-Optic
- Push Button
- CA Burs
- Fits into major implant motors
- 1 Year Warranty

Reg. \$634 **Sale: \$498** Code: SL3019



Maxso® M200 Smart

Non-Optic High Speed Handpiece

- 19W of Torque
- 400,000 RPM
- Push Button Chuck
- Micro-tex Surface
- Triple Water Spray
- 1 Year Warranty

Backend	Code
4 Hole	HP2009

Reg. \$240/ea **Net: \$180/ea**

BUY 3, GET 1 FREE!



Duray® V2

DC X-Ray Unit

- Advanced DC Technology
- Excellent Image Quality
- Built-in DAP allows dose monitoring
- Effortless and stable positioning
- 24 Month Warranty

Code: XR2012

Reg. \$5,250 **Sale: \$4,600**



EZ-Swivel®

Couplers



Beyes® PD 4 Hole
Code: HP6008
Reg. \$152/ea
Net: \$101.33/ea



NSK® Type 4 Hole, Fiber-optic
Code: HP3046
Reg. \$252/ea **Net: \$168/ea**

NSK® Type 6 Hole, Fiber-optic
Code: HP2069
Reg. \$252/ea **Net: \$168/ea**



Kavo® Type 4 Hole, Non-optic
Code: HP2063
Reg. \$252/ea
Net: \$168/ea



Kavo® Type 6 Hole, Fiber-optic
Code: HP2064
Reg. \$252/ea **Net: \$168/ea**

BUY 2, GET 1 FREE!

CareMaster®

Automatic Handpiece Maintenance System

Maintain Your Most Significant Investments in the Office!
Clean 3 Handpieces Per Cycle!

- Maximize the lifespan of handpieces
- Economical, efficient and time-saving
- Filtered air used to rotate the handpieces
- Detachable door for easy cleaning
- Provides care for all brands of handpieces
- 1 Year Warranty

Code: HPL2010
2 Turbines & 1 E-Type Attachment,
Reg. \$1,246

Code: HPL2011
1 Turbines & 2 E-Type Attachments,
Reg. \$1,360

Buy Any Model of CareMaster
Get 1 BeyesLube® 1000ML
& 1 AirLight® M800 Plus FREE



* Requires Airline Connection
* Coupler Sold Separately



Disposable Prophy Angle

BUY 5, GET 1 FREE!

100 Pcs/Box

Soft

Code: 19053

Firm

Code: 19054

Reg. \$46/ea

Net: \$38.33/ea



Maxso® EP12-NS

Hygiene Handpiece

BUY 2, GET 1 FREE!

- 4:1 Nose Cone for DPA
- 1,250–10,000 RPM
- Ceramic Bearings
- Quiet Operation
- Connect Directly to Your Existing Air Motor
- E-Type Connection
- 1 Year Warranty

Code: SL1002

Reg. \$278/ea **Net: \$185.33/ea**



Maxso® Low Speed

E-Type Connection

BUY 3, GET 1 FREE!

- 1:1 Ratio
- 1 Year Warranty

Snap-On

Code: SL2033

Reg. \$137/ea **Net: \$102.75/ea**



Screw-On

Code: SL2034

Reg. \$137/ea **Net: \$102.75/ea**



Latch Type, CA Burs/Latch

Code: SL2028

Reg. \$147/ea **Net: \$110.25/ea**

Push Button, CA Burs/Latch

Code: SL2026

Reg. \$176/ea **Net: \$132/ea**

Straight Nose Cone, HP Burs

Code: ST2010

Reg. \$114/ea **Net: \$85.50/ea**

Air Motor, 4 Holes

Code: MT2010

Reg. \$298/ea **Net: \$223.50/ea**

Mix and Match



Mix & Match any 3. Promotional goods must be of equal or lesser value than the lowest priced item purchased.

GPRO®

Gutta-percha Obturation System

P1 (Heating Pen)

- 0.2 Second Heating Time
- Ergonomic Design
- High Precision Temperature Control
- High-performance Conversion System

Code: GP2003

Reg. \$838- **Sale: \$670**

F1 (Filling Machine)

- Quick Heating within 15 Seconds
- 360° Rotatable Injection Needle
- Ergonomic Design
- Safer Protection Mechanism

Code: GP2002

Reg. \$838- **Sale: \$670**

S1 (Heating Pen + Filling Machine)

Code: GP2001

Reg. \$1,530- **Sale: \$1,225**



ApexPilot® G1

Endo Motor

- Imported Motor 15.5:1 reduction ratio
- Wireless Rechargeable Battery
- Automatic Torque Calibrating
- Nimble Contra-angle

Code: AL2050

Reg. \$690- **Sale: \$555**



ApexPilot® G2

Endo Motor, Connectable to Apex Locator

- Can be connected to matched apex locator Accurator® A2 (sold separately)
- High Performance Brushless Motor
- 6:1 Super Mini Contra-angle
- Efficient versatile functions

Code: AL2060

Reg. \$754- **Sale: \$597**



ApexPilot® G3

Endo Motor with Built-in Apex Locator

- Built-in Apex Locator
- Integrated Length Determination
- All-angle OLED Screen
- Supports most major brands of NiTi files

Code: AL2070

Reg. \$1,258- **Sale: \$999**



Accurator® A2

Apex Locator, Compatible with G2 Endo Motor

- 4.5" wide, color LCD screen provides clear images with high contrast
- Easy-to-use interface with adjustable inclination
- Accurate measurement in wet and dry canals
- Automatic calibration
- Autoclavable file clip, file hook, file probe

Code: AL2080

Reg. \$838- **Sale: \$670**



Comfort-Sonic®

LED Ultrasonic Piezo Scaler

- 3 Modes: Scaling, Endo, Perio
- Built-in Water Regulation
- Daylight LED
- Detachable Handpiece
- Autoclavable Up to 135°C
- Easy Adjustable, Portable Design



P6 Plus

Code: UL2104
Includes: 6 tips
 (EE1, 2 EG1, EG2, EG4, EP1)

Reg. 2 for \$1,224-
Sale: 2 For \$920



P9

Code: UL2103
includes: 8 tips (EE1, EG1, EG4, EG5, EG6, EP1, EP3, EP4D), 1 tip holder & 2 water bottle

Reg. \$1,064-
Sale: \$899

Canaview®

Intraoral Camera

- True Auto-Focus
- Anti-Fog Design
- Responsive Capture Button
- 5 MP High Resolution Images
- Intuitive Camera Lighting System
- No Manual Adjustments Required
- Direct USB Connection to Computer
- Super Light-Weight
- Complimentary Image Management Software

***Compatible with:** Abident, Adstra, Apixia, Apteryx, CaptureLink, Clio Soft, Cuve Hero, DBSWIN, Dentrax Ascend, Dentrax Image, Dexis, Diagora, Eaglesoft, EVAsoft, ImageFX, Kodak, MiPACS, Mogo, Open Dental, Romexis, Schick, Sidexis, Sopro Imaging, Tiger view, Umbie, VaTech EzDent, VisiQuick, Visix, VixWin

Code: EV2060

Reg. \$904/ea
Net: \$452/ea



BUY 1, GET 1 FREE!

BUY 3, GET 1 FREE!

Comfort-Sonic®

Piezo Scaler Handpiece & Accessories

Non-optic:

Compatibility Code
 Beyes & EMS UL2020
Reg. \$227.68/ea Net: \$170.76/ea



Satelec® UL2021
Reg. \$268.09/ea Net: \$201.07/ea

Fiber-optic:

Compatibility Code
 Beyes & EMS UL2022
 Satelec® UL2023



Reg. \$268.09/ea
Net: \$201.07/ea

Tip Wrench:

BUY 3, GET 1 FREE!

WST
Code: UL2420



Reg. \$21.55/ea
Net: \$16.16/ea

E: EMS* & Beyes* Compatible **S:** Satelec* & NSK* Compatible

Scaling Tips:

Used to remove calculus and plaque on supragingival, neck of teeth and hard calculus on the teeth point.

Popular



Used to remove all supragingival calculus and plaque.

Popular



Equivalent to EMS* Type A

Model Code
 E-G1 UL2300
 S-G1 UL2400

Model Code
 E-G4 UL2302
 S-G4 UL2402

Perio Tips:

Used to remove subgingival calculus.

Popular



Used to remove subgingival and deep perio calculus.

Popular



Equivalent to EMS* Type PS

Model Code
 E-P1 UL2305
 S-P1 UL2405

Model Code
 E-P3 UL2306
 S-P3 UL2406

Right angled, to remove calculus from very narrow inter-root spaces and furcation.



Left angled, to remove calculus from very narrow inter-root spaces and furcation.



Model Code
 E-P2R UL2321
 S-P2R UL2431

Model Code
 E-P2L UL2320
 S-P2L UL2430

Reg. \$29.99 - \$42.86/ea

Promo 1:
Buy 5 @ \$25.96/ea

Promo 2:
Buy 10 @ \$21.96/ea

Beyes®

The Next Standard

23-595 Middlefield Rd.,
 Toronto, ON, M1V 3S2, Canada
T: 1-855.603.1888
F: 1-855.720.1228
E: order@beyes.ca
W: www.beyes.ca

Your Dealer for Beyes Products

Scan the QR code for more Beyes products



Printed in Canada
 Document #: ENF045R1
 08.02.23

*Beyes Dental Canada Inc. is not responsible for any typographical errors. Specifications are subject to change without prior notice. All brands are holders of their respective trademarks.